

## Medical Technology

Surgical pads, positioning aids, mattresses, bolus materials and seat cushions made of silicone gel



Hospital/Surgery



Doctor's office



Home care



Domestic care



Baby care



PREMIUM QUALITY



MADE IN GERMANY

Product Catalogue 2024


Sale to commercial customers only. No sale to consumers within the meaning of the § 13 GERMAN CIVIL CODE (BGB)

# Gel-a-Med<sup>®</sup>

Medical sitting and lying



Gel-a-Med<sup>®</sup> products are CE-marked; Regulation (EU) 2017/745 (MDR) and belong to medical device class I. Gel-a-Med<sup>®</sup> Certification according to 13485:2021, 9001:2015, 14001:2015



## Product features of Gel-a-Med® seat cushions, surgical pads and positioning aids

Made in Germany. For more than 30 years, we have been manufacturing silicone gel products for you to the highest medical standards at our site in Düsseldorf. Gel-a-Med® products are used in hospitals, children's hospitals, outpatient surgery centers, rehabilitation clinics, senior citizen facilities, physiotherapy practices and in home care. Compliance with quality and environmental management in accordance with ISO 9001:2015, ISO 14001:2015 and ISO 13485, the internationally relevant standard for quality management systems of medical device manufacturers, as well as the European guidelines according to MDR guarantee first-class product quality.

### Valuable properties

- ✓ High elasticity
- ✓ Optimum adaptation to body contours
- ✓ Even distribution of body pressure
- ✓ Cushioning, absorbing effect
- ✓ No sink in of body for maximum elbow room
- ✓ Form stability, temperature resistance
- ✓ Long-term flexibility
- ✓ Long shelf life
- ✓ Adaptation to body temperature
- ✓ Easy to clean and sanitize

### Possible applications

- ✓ For decubitus prophylaxis and therapy
- ✓ For efficient pressure relief
- ✓ For prevention of circulatory failure and nerve damage
- ✓ Highly convenient and comfortable positioning of sensitive body parts
- ✓ For general well-being
- ✓ For long time sitting and lying

### The medical plus

- ✓ Made of high quality, health uncritical silicone
- ✓ Good skin compatibility
- ✓ Free of allergens, latex, harmful substances or plasticizers
- ✓ Leak-proof due to special molecule composition
- ✓ Up to 90 percent X-ray permeability

Seating cushions and toilet cushions  
see page 17-18



IF APPLICABLE IN YOUR COUNTRY!

### Apply for aids from statutory health insurance!

Gel-a-Med® seat cushions and toilet cushions can be applied for as aids from statutory health insurance


#### All-silicone products

Aids-Number:

11.39.02.0001

## Gel-a-Med® pads ensure optimal patient positioning in the OR


During surgery, the body is exposed to particularly high pressure loads. Pressure ulcers as a result of prolonged positioning or shear forces develop deep within the tissue and often become visible only days later. Pressure ulcer prophylaxis during surgery through proper positioning and effective pressure relief is therefore necessary from a medical, nursing and economic point of view.





## Gel-a-Med® surgical pads and positioning aids offer ...

- ✓ Best possible access to the surgical area
- ✓ Prevention of pressure sores/decubitus during and after surgery
- ✓ Prevention of circulatory disorders
- ✓ Prevention of nerve damage
- ✓ Easy handling



## Advantages of Gel-a-Med® products made of pressure relieving silicone gel



### Quality from the leading manufacturer in Germany

Gel-a-Med® medical products made of high-quality silicone gel have been produced entirely in Germany for 30 years - for clinics, outpatient surgical centers, nursing homes and the home care sector. Trust in the experience of the German market leader.



### Highly flexible

Due to the unique molecular structure of the silicone gel, Gel-a-Med® medical products are highly elastic, soft and pliable and they reduce the shearing forces acting on the body to a minimum. Furthermore, they permanently retain their shape and they are leak-proof.



### Certified medical product

Gel-a-Med® medical products are manufactured in compliance with high legal standards. Gel-a-Med GmbH is certified according to DIN ISO 13485 for Class I medical products.



### Skin compatible and antiallergic

Gel-a-Med® medical products are characterized by their good skin compatibility. They are free of allergens, plasticizers, latex and harmful substances.




### Suitable for positioning sensitive body parts

Gel-a-Med® medical products offer maximum protection when lying and sitting as well as during transport. The cushioning effect enables safe and comfortable positioning. An enlarged contact surface reduces and absorbs local pressure - the skin is being unburdened.



### Breathable

Gel-a-Med® medical products are encased in a breathable membrane. Microscopic holes in the sturdy membrane provide a comfortable and medically favorable microclimate.



PREMIUM QUALITY  
MADE IN GERMANY



## Gel-a-Med medical products for an effective and efficient decubitus prophylaxis

Any limitation of movement leads to an increased pressure load, which causes local damage to the skin and the underlying tissue.

Gel-a-Med® medical products are indispensable aids for positioning patients at risk of decubitus in hospitals, outpatient surgery centers, nursing homes, rehabilitation centers and in the home care sector.



### Temperature resistant

Gel-a-Med® medical products are temperature resistant from -20 °C to +55 °C and adapt optimally to body temperature as well as store heat. Gel-a-Med® products remain permanently elastic and pliable even at low temperatures.



### Resistant to liquids, water and steam

Gel-a-Med® medical products are liquid repellent, steam resistant and quick to dry after use.



### Disinfectable

Gel-a-Med® products can be easily sanitized with disinfectants. Further information can be found in the instructions for use at [www.gel-a-med.com](http://www.gel-a-med.com)



### Entirely odorless

Gel-a-Med® medical products remain odorless when used as intended.




### Sustainable

Due to their long-lasting elasticity and structural stability as well as their good sanitizability, Gel-a-Med® medical products can be used over a long period of time, thus conserving resources.



# Surgical pads




Operating table pads  
in various sizes

## Gel-a-Med® Surgical pads · Operating table pads

Silicone gel

Article-No.	Description	Dimensions W x D x H	Weight
M1-1608008	Surgical pad	16 x 8 x 0,8 cm	0,10 kg
M1-1711015	Surgical pad	17 x 11 x 1,5 cm	0,25 kg
M1-2020010	Surgical pad	20 x 20 x 1,0 cm	0,40 kg
M1-3008008	Surgical pad	30 x 8 x 0,8 cm	0,20 kg
M1-3020010	Surgical pad	30 x 20 x 1,0 cm	0,60 kg
M1-3030010	Surgical pad	30 x 30 x 1,0 cm	0,90 kg
M1-3030020	Surgical pad	30 x 30 x 2,0 cm	1,80 kg
M1-4012010	Surgical pad	40 x 12 x 1,0 cm	0,45 kg
M1-4020018	Surgical pad	40 x 20 x 1,8 cm	1,40 kg
M1-4025010	Surgical pad	40 x 25 x 1,0 cm	1,00 kg
M1-4040010	Surgical pad	40 x 40 x 1,0 cm	1,60 kg
M1-4040020	Surgical pad	40 x 40 x 2,0 cm	3,00 kg
M1-4342020	Surgical pad	43 x 42 x 2,0 cm	3,60 kg
M1-4542020	Surgical pad	45 x 42 x 2,0 cm	3,80 kg
M1-4842020	Surgical pad	48 x 42 x 2,0 cm	4,00 kg
M1-5020010	Surgical pad	50 x 20 x 1,0 cm	1,00 kg
M1-5230010	Surgical pad	52 x 30 x 1,0 cm	1,50 kg
M1-5252010	Surgical pad	52 x 52 x 1,0 cm	2,50 kg
M1-5252020	Surgical pad	52 x 52 x 2,0 cm	5,40 kg
M1-6012010	Surgical pad	60 x 12 x 1,0 cm	0,70 kg



Article-No.	Description	Dimensions W x D x H	Weight
M1-6020010	Surgical pad	60 x 20 x 1,0 cm	1,20 kg
M1-6040010	Surgical pad	60 x 40 x 1,0 cm	2,40 kg
M1-6040015	Surgical pad	60 x 40 x 1,5 cm	3,60 kg
M1-6040020	Surgical pad	60 x 40 x 2,0 cm	4,60 kg
M1-7550013	Surgical pad	75 x 50 x 1,3 cm	4,60 kg
M1-8060015	Surgical pad	80 x 60 x 1,5 cm	7,20 kg
M1-9050015	Surgical pad	90 x 50 x 1,5 cm	6,80 kg
M1-9050020	Surgical pad	90 x 50 x 2,0 cm	9,00 kg
M1-9060015	Surgical pad	90 x 60 x 1,5 cm	8,10 kg
M1-10050010	Surgical pad	100 x 50 x 1,0 cm	5,00 kg
M1-11050015	Surgical pad	110 x 50 x 1,5 cm	8,20 kg
M1-11050020	Surgical pad	110 x 50 x 2,0 cm	11,00 kg
M1-11060020	Surgical pad	110 x 60 x 2,0 cm	13,20 kg
M1-11552010	Surgical pad	115 x 52 x 1,0 cm	6,00 kg
M1-12070020	Surgical pad	120 x 70 x 2,0 cm	16,80 kg
M1-17550015	Surgical pad	175 x 50 x 1,5 cm	13,00 kg
M1-18052010	Surgical pad	180 x 52 x 1,0 cm	9,40 kg
M1-19050015	Surgical pad	190 x 50 x 1,5 cm	14,20 kg
M1-5252010V	Surgical pad with sacral cutout	52 x 52 x 1,0 cm	2,50 kg
M1-9050010V	Surgical pad with sacral cutout	90 x 50 x 1,0 cm	4,50 kg




## Special surgical pads



### Gel-a-Med® Perforated surgical pads

Silicone gel

Perforated for good penetration of water and water vapor, for treatment wards for severe burn injuries, pathology or morticians



Perforated surgical pad

Special color

Article-No.	Description	Dimensions W x D x H	Weight
M1-45090020P	Perforated surgical pads	45 x 90 x 2,0 cm	6,70 kg
M1-19052020P	Perforated surgical pads	190 x 52 x 2,0 cm	17,50 kg



# Cold and heat therapy pads



Thermotherapy pad  
in various sizes

## Gel-a-Med® Cold and heat therapy pads

Silicone gel

Temperature resistant from -20 °C to +55 °C, heat retentive

Article-No.	Description	Dimensions W x D x H	Weight
M1-1608008T	Thermotherapy pad	16 x 8 x 0,8 cm	0,10 kg
M1-1711015T	Thermotherapy pad	17 x 11 x 1,5 cm	0,25 kg
M1-2020010T	Thermotherapy pad	20 x 20 x 1,0 cm	0,40 kg
M1-3020010T	Thermotherapy pad	30 x 20 x 1,0 cm	0,60 kg
M1-3030010T	Thermotherapy pad	30 x 30 x 1,0 cm	0,90 kg
M1-3030020T	Thermotherapy pad	30 x 30 x 2,0 cm	1,80 kg
M1-4012010T	Thermotherapy pad	40 x 12 x 1,0 cm	0,45 kg
M1-4020018T	Thermotherapy pad	40 x 20 x 1,8 cm	1,40 kg
M1-4025010T	Thermotherapy pad	40 x 25 x 1,0 cm	1,00 kg
M1-4040010T	Thermotherapy pad	40 x 40 x 1,0 cm	1,60 kg
M1-4040020T	Thermotherapy pad	40 x 40 x 2,0 cm	3,00 kg
M1-4342020T	Thermotherapy pad	43 x 42 x 2,0 cm	3,60 kg
M1-4542020T	Thermotherapy pad	45 x 42 x 2,0 cm	3,80 kg
M1-4842020T	Thermotherapy pad	48 x 42 x 2,0 cm	4,00 kg
M1-5020010T	Thermotherapy pad	50 x 20 x 1,0 cm	1,00 kg
M1-5230010T	Thermotherapy pad	52 x 30 x 1,0 cm	1,50 kg
M1-5252010T	Thermotherapy pad	52 x 52 x 1,0 cm	2,50 kg
M1-5252020T	Thermotherapy pad	52 x 52 x 2,0 cm	5,40 kg
M1-6012010T	Thermotherapy pad	60 x 12 x 1,0 cm	0,70 kg
M1-6020010T	Thermotherapy pad	60 x 20 x 1,0 cm	1,20 kg
M1-6040010T	Thermotherapy pad	60 x 40 x 1,0 cm	2,40 kg
M1-6040015T	Thermotherapy pad	60 x 40 x 1,5 cm	3,60 kg
M1-6040020T	Thermotherapy pad	60 x 40 x 2,0 cm	4,60 kg




## Positioning head / face



### Gel-a-Med® Forehead cushion · Face cushion

Silicone gel



Forehead cushion / Face cushion

Article-No.	Description	Dimensions W x D x H	Weight
M1-1608008S	Forehead cushion / Face cushion	16 x 8 x 0,8 cm	0,10 kg
M1-3008008S	Forehead cushion / Face cushion	30 x 8 x 0,8 cm	0,20 kg



### Gel-a-Med® Head rings · Head cushion


Silicone gel



Head ring, closed



Head ring, open



Head cushion large, closed

Article-No.	Description	Dimensions W x D x H	Weight
M1-1313040K	Head ring small, closed Ø 6 cm	13 x 13 x 4,0 cm	0,31 kg
M1-1313040K-1	Head ring small, open Ø 6 cm	13 x 13 x 4,0 cm	0,27 kg
M1-2121040K	Head ring large, closed Ø 6 cm	21 x 21 x 4,0 cm	1,30 kg
M1-2121040K-1	Head ring large, open Ø 6 cm	21 x 21 x 4,0 cm	1,20 kg
M1-3030040K	Head cushion large, closed Ø 6 cm	30 x 30 x 4,0 cm	2,70 kg



### Gel-a-Med® Head cushion anatomical prone position

1,5 cm silicone gel, foam and PU cover

For positioning the head in the prone position in the OR



Head cushion anatomical

Article-No.	Description	Cut-out	Weight
M1-2621120KA	Head cushion anatomical	18 x 12 cm	0,60 kg

# Positioning head / face



## Gel-a-Med® Head rings round

For positioning the head in the operating or nursing bed



Head ring round, closed



Head ring, open half moon

Article-No.	Description	Outer Ø	Inner Ø	Strength	Weight
M1-1706025K	Head ring round, closed	17,0 cm	6,0 cm	2,5 cm	0,35 kg
M1-205045030K	Head ring round, closed	20,5 cm	4,5 cm	3,0 cm	0,62 kg
M1-20510020K	Head ring round, closed	20,5 cm	10,0 cm	2,0 cm	0,37 kg
M1-19060025KH	Head ring, open half moon	19,0 cm	6,0 cm	2,5 cm	0,30 kg



## Gel-a-Med® Head rings round with flat contact surface

Color blue, flat support surface, to protect the patient's head and ears in the lateral position and easy access for intubation



Head ring, open half moon

Article-No.	Description	Outer Ø	Inner Ø	Strength	Weight
M1-200950405KOP	Head ring, open half moon	20,0 cm	9,5 cm	4,5 cm	0,55 kg



## Gel-a-Med® Head rings round with flat contact surface

Color blue, flat support surface, for positioning the head in the operating or nursing bed in supine position



Head ring round



Head ring large, closed with dent

Article-No.	Description	Outer Ø	Inner Ø	Strength	Weight
M1-20545050KPD	Head ring large, closed with dent	20,5 cm	4,5 cm	5,0 cm	1,35 kg
M1-200950405KP	Head ring round	20,0 cm	9,5 cm	4,5 cm	0,85 kg




## Positioning arm / elbow



### Gel-a-Med® Arm / Elbow positioning

Silicone gel

Protection and support for the sensitive arm-elbow area



Article-No.	Description	Dimensions W x D x H	Weight
M1-1711015A	Arm / Hand positioning pad	17 x 11 x 1,5 cm	0,25 kg
M1-3020010A	Arm positioning pad	30 x 20 x 1,0 cm	0,60 kg
M1-4012010A	Arm positioning pad	40 x 12 x 1,0 cm	0,45 kg
M1-4025010A	Arm positioning pad	40 x 25 x 1,0 cm	1,00 kg
M1-5020010A	Arm positioning pad	50 x 20 x 1,0 cm	1,00 kg
M1-5020010A-K	Arm positioning pad w/ velcro strap	50 x 20 x 1,0 cm	1,00 kg
M1-6012010A	Arm positioning pad	60 x 12 x 1,0 cm	0,70 kg
M1-6020010A	Arm positioning pad	60 x 20 x 1,0 cm	1,20 kg
M1-5115035AK	Arm positioning pad concavev	51 x 15 x 3,5 cm	2,30 kg
M1-5115060AK	Arm positioning pad concavev w/ velcro strap	51 x 15 x 6,0 cm	4,59 kg
M1-6115050AK	Arm positioning pad concavev w/ velcro strap	61 x 15 x 5,0 cm	4,58 kg

# Positioning leg / knee / heel



## Gel-a-Med® Leg positioning pads · Leg shell-pads

Silicone gel



Leg positioning pad

Article-No.	Description	Dimensions W x D x H	Weight
M1-4025010B	Leg positioning pad	40 x 25 x 1,0 cm	1,00 kg
M1-4020018B	Leg positioning pad	40 x 20 x 1,8 cm	1,40 kg
M1-5230010B	Leg positioning pad	52 x 30 x 1,0 cm	1,50 kg
M1-6020010B	Leg positioning pad	60 x 20 x 1,0 cm	1,20 kg



## Gel-a-Med® Knee pad Silicone gel Knee pad half roll Foam core RG 33/60, silicone gel topper



Knee pad silicone gel

Knee pad half roll  
foam core

Article-No.	Description	Dimensions W x D x H	Weight
M1-5115050KP	Knee pad, silicone gel	51 x 15 x 5,0 cm	3,00 kg
M3-5015100KP	Knee pad half roll, silicone/foam	50 x 15 x 10,0 cm	1,20 kg



## Gel-a-Med® Heel pad

Silicone gel



Heel pad, in pairs

Article-No.	Description	Dimensions W x D x H	Weight
M1-1508040F	Heel pad, in pair	15 x 8 x 4,0 cm	0,80 kg
M1-2010050F	Heel pad, in pair	20 x 10 x 5,0 cm	1,60 kg



## Positioning leg / heel



### Gel-a-Med® Concave leg protector

Serves to securely support and position of the leg and heel area

Article-No.	Description	Dimensions W x D x H	Weight
M1-52180400BK	Concave leg protector	52 x 18 x 4,0 cm	2,85 kg



### Gel-a-Med® Leg protector

Serves to securely support and position of the leg area

Article-No.	Description	Dimensions W x D x H	Weight
M1-38210425B	Leg protector	38 x 21 x 2,5 cm	2,00 kg



# Positioning



## Gel-a-Med® Positioning aid half-roll

Foam core RG 33/60, silicone gel topper



Positioning aids half-roll in various sizes

Article-No.	Description	Dimensions W x D x H	Weight
M3-3015100HR	Positioning aid half-roll	30 x 15 x 10 cm	0,70 kg
M3-3015150HR	Positioning aid half-roll	30 x 15 x 15 cm	0,90 kg
M3-4015100HR	Positioning aid half-roll	40 x 15 x 10 cm	0,80 kg
M3-4015150HR	Positioning aid half-roll	40 x 15 x 15 cm	1,20 kg
M3-5015100HR	Positioning aid half-roll	50 x 15 x 10 cm	1,00 kg
M3-5015150HR	Positioning aid half-roll	50 x 15 x 15 cm	1,50 kg



## Gel-a-Med® Positioning aid / Positioning cushion

Silicone gel



Positioning cushions in various sizes

Article-No.	Description	Dimensions W x D x H	Weight
M1-3010050	Positioning cushion	30 x 10 x 5,0 cm	1,20 kg
M1-4010050	Positioning cushion	40 x 10 x 5,0 cm	1,60 kg
M1-4020040	Positioning cushion	40 x 20 x 4,0 cm	3,00 kg
M1-5010050	Positioning cushion	50 x 10 x 5,0 cm	2,00 kg
M1-5115050	Positioning cushion	51 x 15 x 5,0 cm	3,00 kg



## Gel-a-Med® Side support cushion

With 1 cm silicone gel pad

For stabilizing and positioning individual parts of the body in a lateral position



Side support cushion

Article-No.	Description	Dimensions W x D x H	Weight
M3-7250160SL	Side support cushion	72 x 50 x 16 cm	3,75 kg



## Positioning



### Gel-a-Med® Belly pad large

With 1 cm silicone gel pad

For support during operations in the prone position



Belly pad large

Article-No.	Description	Dimensions W x D x H	Weight
M3-11850080BP	Belly pad large	118 x 50 x 16 cm	5,90 kg



### Gel-a-Med® Belly pad with interchangeable wedge

With 1 cm silicone gel pad

For support during operations in the prone position with interchangeable wedge and Velcro



Belly pad with interchangeable wedge

Article-No.	Description	Dimensions W x D x H	Weight
M3-7550080BP	Belly pad with interchangeable wedge	75 x 50 x 8 cm	3,75 kg



### Gel-a-Med® Renal puncture cushion with belt

With 1 cm silicone gel pad



Renal puncture cushion with belt

Article-No.	Description	Dimensions W x D x H	Weight
M3-8154021NP	Renal puncture cushion with belt	81 x 54 x 21 cm	4,00 kg

# Positioning



## Gel-a-Med® Positioning cushion

Including white cotton cover and 2 carrying handles, for stabilizing and positioning individual parts of the body



Positioning cushion

with cotton cover


Article-No.	Description	Dimensions W x D x H	Weight
M1-4343120	Positioning cushion	43 x 43 x 12 cm	22,00 kg



## Gel-a-Med® Positioning roller

Silicone gel

For stabilizing and positioning individual parts of the body



Positioning roller

Article-No.	Description	Ø	Length	Weight
M3-5015000LR	Positioning roller	15 cm	50 cm	3,00 kg



## Gel-a-Med® Half roll

Blue silicone gel

For stabilizing and positioning individual parts of the body



Half roll

Article-No.	Description	Dimensions W x D x H	Weight
M3-50130090HR	Half roll	50 x 13 x 9 cm	4,00 kg




## Positioning



### Gel-a-Med® Positioning roller

With 1 cm silicone gel insert and foam core

Including white PU cover, for stabilizing and positioning individual parts of the body



Positioning roller

Article-No.	Description	Ø	Length	Weight
M3-4022000LR-PU	Positioning roller	22 cm	40 cm	3,00 kg



### Gel-a-Med® Surgical positioning roll

Foam core without silicone gel

With white PU cover, for stabilizing and positioning individual parts of the body



Surgical positioning roll with cover

Article-No.	Description	Ø	Length	Weight
M3-4020LRS-PU	Surgical positioning roll	20 cm	40 cm	1,20 kg

# Positioning



## Gel-a-Med® Positioning cube

Foam

To relieve the spine and decongest the legs, PU cover washable at 95 degrees, resistant to water, blood and urine and suitable for wipe disinfection



Positioning cube

Article-No.	Description	Dimensions W x D x H	Weight
M3-303030SW-B	Positioning cube blue	30 x 30 x 30 cm	3,50 kg
M3-303030SW-G	Positioning cube gray	30 x 30 x 30 cm	3,50 kg
M3-303030SW-Y	Positioning cube yellow	30 x 30 x 30 cm	3,50 kg
M3-303030SW-S	Positioning cube black	30 x 30 x 30 cm	3,50 kg
M3-303030SW-W	Positioning cube white	30 x 30 x 30 cm	3,50 kg
M3-404040SW-B	Positioning cube blue	40 x 40 x 40 cm	4,00 kg
M3-404040SW-G	Positioning cube gray	40 x 40 x 40 cm	4,00 kg
M3-404040SW-Y	Positioning cube yellow	40 x 40 x 40 cm	4,00 kg
M3-404040SW-S	Positioning cube black	40 x 40 x 40 cm	4,00 kg
M3-404040SW-W	Positioning cube white	40 x 40 x 40 cm	4,00 kg
M3-505050SW-B	Positioning cube blue	50 x 50 x 50 cm	4,50 kg
M3-505050SW-G	Positioning cube gray	50 x 50 x 50 cm	4,50 kg
M3-505050SW-Y	Positioning cube yellow	50 x 50 x 50 cm	4,50 kg
M3-505050SW-S	Positioning cube black	50 x 50 x 50 cm	4,50 kg
M3-505050SW-W	Positioning cube white	50 x 50 x 50 cm	4,50 kg



## Positioning



### Gel-a-Med® Stepped positioning cube (3 steps)

Foam

To relieve the spine and decongest the legs, PU cover washable at 95 degrees, resistant to water, blood and urine and suitable for wipe disinfection.



Stepped positioning cube

Article-No.	Description	Dimensions W x D x H	Weight
M3-504540LW-B	Stepped positioning cube blue	50 x 45 x 40 cm	4,00 kg
M3-504540LW-G	Stepped positioning cube gray	50 x 45 x 40 cm	4,00 kg
M3-504540LW-Y	Stepped positioning cube yellow	50 x 45 x 40 cm	4,00 kg
M3-504540LW-S	Stepped positioning cube black	50 x 45 x 40 cm	4,00 kg
M3-504540LW-W	Stepped positioning cube white	50 x 45 x 40 cm	4,00 kg



### Gel-a-Med® Stepped positioning cube adjustable (2-piece)

Foam

to relieve the spine and decongest the legs, PU cover washable at 95 degrees, resistant to water, blood and urine and suitable for wipe disinfection. Adjustable in height from 50-75 cm using Velcro tape, 2-piece.




Stepped positioning cube

Article-No.	Description	Dimensions W x D x H	Weight
M3-505075LWV-B	Stepped positioning cube blue	50 x 50 x 50-75 cm	4,60 kg
M3-505075LWV-G	Stepped positioning cube gray	50 x 50 x 50-75 cm	4,60 kg
M3-505075LWV-Y	Stepped positioning cube yellow	50 x 50 x 50-75 cm	4,60 kg
M3-505075LWV-S	Stepped positioning cube black	50 x 50 x 50-75 cm	4,60 kg
M3-505075LWV-W	Stepped positioning cube white	50 x 50 x 50-75 cm	4,60 kg



# Bed toppers / Mattresses



## Gel-a-Med® Bed toppers

Silicone gel


Including cotton linen, white

Article-No.	Description	Dimensions W x D x H	Weight
M1-12070013B	Bed topper	120 x 70 x 1,3 cm	10,00 kg
M1-20090010B	Bed topper	200 x 90 x 1,0 cm	18,00 kg
M1-20090013B	Bed topper	200 x 90 x 1,3 cm	22,20 kg
M2-6040013B	Bed topper	60 x 40 x 1,3 cm	3,12 kg



## Gel-a-Med® Mattresses

Healthcare/ nursing mattresses with gel inlay  
PU cover water repellent, removable, washable  
at 60 degrees and suitable for dryer



Article-No.	Description	Dimensions W x D x H	Weight
M6-90200140	Healthcare mattress	90 x 200 x 14 cm	15,00 kg
M6-100200140	Healthcare mattress	100 x 200 x 14 cm	18,00 kg



## Gel-a-Med® Evacuation mattress

4,0 cm foam + 0,8 cm silicone gel layer

8 carrying handles with silicone gel layer (4 on each side), 2 carrying straps, 3 fastening straps with velcro, for patient transport. PU cover yellow, disinfectable, blood- urine- and water-resistant. Can be easily folded for transport, carried like a suitcase and fixed with Velcro straps.



Article-No.	Description	Dimensions W x D x H	Weight
M1-19555050EM	Evacuation mattress	195 x 55 x 5 cm	9,60 kg



## Radiation therapy



### Gel-a-Med® Bolus mat / Positioning aid radiographic / Radiotherapy mats

Transparent silicone gel

Suitable for radio/radiation therapy



Rectangular bolus

Article-No.	Description	Dimensions W x D x H	Weight
M1-1010005ST-T	Bolus mat XS, squared	10 x 10 x 0,5 cm	0,05 kg
M1-1520005ST-T	Bolus mat S, rectangular	15 x 20 x 0,5 cm	0,15 kg
M1-3030005ST-T	Bolus mat M, squared	30 x 30 x 0,5 cm	0,45 kg
M1-2040005ST-T	Bolus mat M, rectangular	20 x 40 x 0,5 cm	0,40 kg
M1-4040005ST-T	Bolus mat L, squared	40 x 40 x 0,5 cm	0,80 kg
M1-5050005ST-T	Bolus mat XL, squared	50 x 50 x 0,5 cm	1,25 kg



### Gel-a-Med® Bolus cushion / Breast positioning cushion

Blue silicone gel

Suitable for radio/radiation therapy



Breast positioning cushion

Article-No.	Description	Dimensions W x D x H	Weight
M3-50130090BK	Breast positioning cushion	50 x 13 x 9 cm	5,85 kg

# Seat cushion / Wheelchair cushion



Seat cushion height 2,0 cm

Seat cushion height 4,0 cm

Seat cushion with linen

**Gel-a-Med®**

## Seat cushion / Wheelchair cushion


Silicone gel

Including quilted cotton linen black, with handle, Oeko-Tex® Standard 100

Number for aids from statutory health insurance 11.39.02.0001

Article-No.	Description	For body weight	Dimensions W x D x H	Weight
M2-4040020	Wheelchair cushion	up to 70 kg	40 x 40 x 2,0 cm	3,00 kg
M2-4342020	Wheelchair cushion	up to 70 kg	43 x 42 x 2,0 cm	3,60 kg
M2-4542020	Wheelchair cushion	up to 70 kg	45 x 42 x 2,0 cm	3,80 kg
M2-4842020	Wheelchair cushion	up to 70 kg	48 x 42 x 2,0 cm	4,00 kg
M2-4040025	Wheelchair cushion	up to 90 kg	40 x 40 x 2,5 cm	3,90 kg
M2-4342025	Wheelchair cushion	up to 90 kg	43 x 42 x 2,5 cm	4,40 kg
M2-4542025	Wheelchair cushion	up to 90 kg	45 x 42 x 2,5 cm	4,60 kg
M2-4842025	Wheelchair cushion	up to 90 kg	48 x 42 x 2,5 cm	4,90 kg
M2-4040035	Wheelchair cushion	up to 120 kg	40 x 40 x 3,5 cm	5,60 kg
M2-4342035	Wheelchair cushion	up to 120 kg	43 x 42 x 3,5 cm	6,00 kg
M2-4542035	Wheelchair cushion	up to 120 kg	45 x 42 x 3,5 cm	6,60 kg
M2-4842035	Wheelchair cushion	up to 120 kg	48 x 42 x 3,5 cm	6,80 kg
M2-5050035	Wheelchair cushion	up to 120 kg	50 x 50 x 3,5 cm	8,40 kg
M2-4040040	Wheelchair cushion	from 120 kg	40 x 40 x 4,0 cm	6,10 kg
M2-4342040	Wheelchair cushion	from 120 kg	43 x 42 x 4,0 cm	6,90 kg
M2-4542040	Wheelchair cushion	from 120 kg	45 x 42 x 4,0 cm	7,20 kg
M2-4842040	Wheelchair cushion	from 120 kg	48 x 42 x 4,0 cm	7,70 kg
M2-5050040	Wheelchair cushion	from 120 kg	50 x 50 x 4,0 cm	9,50 kg

## Which seat cushion thickness for which body weight?





## Seat cushion



### Gel-a-Med® Seat wedge cushion

Silicone gel

Wedge-shaped, for stabilizing and positioning individual parts of the body. Cover included.  
Number for aids from statutory health insurance 11.39.02.0001



Seat wedge cushion without cover

Article-No.	Description	Dimensions W x D x H	Weight
M1-3242030SK	Seat wedge cushion without cover	32 x 42 x 3,0 cm	2,40 kg
M2-4040040SK	Seat wedge cushion	40 x 40 x 4,0 cm	3,60 kg
M2-4342040SK	Seat wedge cushion	43 x 42 x 4,0 cm	4,10 kg
M2-4542040SK	Seat wedge cushion	45 x 42 x 4,0 cm	4,30 kg
M2-4842040SK	Seat wedge cushion	48 x 42 x 4,0 cm	4,60 kg



### Gel-a-Med® Seat cushion/Ring cushion

Silicone gel

Including PU cover gray, up to body weight 200 kg, used for local therapy and prophylaxis for pressure sores, pain in the buttock region, after operations in the perineal area, for anal lesions, hemorrhoids or during or after pregnancy | Number for aids from statutory health insurance 11.39.02.0001



Seat cushion/Ring cushion



with cover

Article-No.	Description	Outer Ø	Inner Ø	Strength	Weight
M2-3838060R-PU	Seat cushion/Ring cushion	38 cm	12,5 cm	6,0 cm	5,00 kg
M2-4040006R-PU	Seat cushion/Ring cushion	40 cm	12,5 cm	6,0 cm	5,80 kg
M2-4343060R-PU	Seat cushion/Ring cushion	43 cm	12,5 cm	6,0 cm	6,80 kg
M2-4545060R-PU	Seat cushion/Ring cushion	45 cm	12,5 cm	6,0 cm	7,40 kg

# Toilet cushion / Shower seat cushion




Gel-a-Med®

Toilet seat cushion · Shower seat cushion, closed

Silicone gel

Including quilted cotton linen black, with handle, Oeko-Tex® Standard 100

For body weight up to 70 kg | Number for aids from statutory health insurance 11.39.02.0001



Toilet cushion, closed

with linen

anatomical, closed

Article-No.	Description	Inner Ø	Dimensions W x D x H	Weight
M2-4040020T	Toilet cushion, closed	18 cm	40 x 40 x 2,0 cm	2,10 kg
M2-4342020T	Toilet cushion, closed	18 cm	43 x 42 x 2,0 cm	2,50 kg
M2-4542020T	Toilet cushion, closed	18 cm	45 x 42 x 2,0 cm	2,70 kg
M2-4342025TA-PU	Toilet cushion anatomical, closed	18 x 25 cm	43 x 42 x 2,5 cm	2,70 kg
M2-4040025TA-PU	Toilet cushion anatomical, closed	18 x 25 cm	40 x 40 x 2,5 cm	2,40 kg




Gel-a-Med®

Toilet cushion · Shower seat cushion, open

Silicone gel

Including quilted cotton linen black, with handle, Oeko-Tex® Standard 100

For body weight up to 70 kg | Number for aids from statutory health insurance 11.39.02.0001



Toilet cushion, open

with linen

anatomical, open (8 cm)

Article-No.	Description	Inner Ø	Dimensions W x D x H	Weight
M2-4040020T-1	Toilet cushion, open	18 cm	40 x 40 x 2,0 cm	1,90 kg
M2-4342020T-1	Toilet cushion, open	18 cm	43 x 42 x 2,0 cm	2,20 kg
M2-4542020T-1	Toilet cushion, open	18 cm	45 x 42 x 2,0 cm	2,40 kg
M2-4342025TA-PU-1	Toilet cushion anatomical, open	18 x 25 cm	43 x 42 x 2,5 cm	2,70 kg
M2-4040025TA-PU-1	Toilet cushion anatomical, open	18 x 25 cm	40 x 40 x 2,5 cm	2,30 kg



Gel-a-Med®

Universal toilet cushion for standard toilets

Silicone gel

Fits on all common standard toilets, oval, anti-slip underside. Application for therapy and

prophylaxis of decubitus ulcers | Number for aids from statutory health insurance 11.39.02.0001



Universal toilet cushion

Article-No.	Description	Cut-out	Dimensions W x D x H	Weight
M1-3945020TG	Universal toilet cushion	19 x 27 cm	39 x 45 x 2,0 cm	1,80 kg



## Infants / Babies



### Gel-a-Med® Positioning cushion, infants

For stabilizing and positioning individual parts of the body in the incubator, operating or nursing bed



Positioning cushion round S

Article-No.	Description	Ø	Strength	Weight
M1-1000015LKS	Positioning cushion round S, infants	10 cm	1,5 cm	0,08 kg



### Gel-a-Med® Positioning cushion, infants

For stabilizing and positioning individual parts of the body in the incubator, operating or nursing bed



Positioning cushion 8 angular



Positioning cushion hexagon M



Positioning cushion bone M

Article-No.	Description	Dimensions W x D x H	Weight
M1-1010015LK	Positioning cushion 8 angular, infants	10 x 10 x 1,5 cm	0,08 kg
M1-1918020LKM	Positioning cushion hexagon M, infants from 2,5 kg	19 x 18 x 2,0 cm	0,40 kg
M1-2215015LKM	Positioning cushion bone M, infants from 2,5 kg	22 x 12-15 x 1,5 cm	0,30 kg



### Gel-a-Med® Head ring round, closed, infants

For positioning the head in the incubator, operating or nursing bed



Head ring round, closed

Article-No.	Description	Dimensions W x D x H	Weight
M1-0904015K	Head ring round, closed, infants	9 x 4 x 1,5 cm	0,05 kg



## Infants / Babies



### Gel-a-Med® Changing table positioning aid Light

With 2 cm PUR foam and 1 cm silicone gel

Side sections 70 x 14 x 3 cm and head section 44 x 14 x 3 cm made of PUR foam, foldable, for positioning individual body parts, for neonatal and children's ward. Including PU cover with zipper.



Changing table positioning aid Light

Article-No.	Description	Lying area W x D x H	Weight
M3-7044030WL	Changing table positioning aid Light	70 x 44 x 3,0 cm	2,80 kg



### Gel-a-Med® Changing table mat

Silicone gel

For positioning individual body parts, for neonatal and children's ward. Including PU cover with zipper



Changing table mat

Article-No.	Description	Dimensions W x D x H	Weight
M2-10060010W	Changing table mat	100 x 60 x 1,0 cm	6,00 kg



### Gel-a-Med® Warming bed pad

Silicone gel

Temperature-resistant from -20 °C to +55 °C, heat-retaining, for positioning individual body parts in incubators, operating or nursing beds



Warming bed pad

Article-No.	Description	Dimensions W x D x H	Weight
M1-4974013W	Warming bed pad	49 x 74 x 1,3 cm	4,70 kg



## Supply blankets



### Gel-a-Med® Supply blanket

100 % polyester fleece thermally bonded

Flame retardant according to BS 5852, Oeko-Tex Standard 100, antibacterial, skin-friendly, color white, PU 10 pieces (individually packed in PE bag). For use in the RTW or for covering patients before and after surgical procedures.



Supply blanket

Article-No.	Description	Dimensions W x D x H	Weight	PU
M7-110198010V	Supply blanket	110 x 198 x 1,0 cm	2,20 kg	10 pieces
M7-130198010V	Supply blanket	130 x 198 x 1,0 cm	2,40 kg	10 pieces
M7-150198010V	Supply blanket	150 x 198 x 1,0 cm	2,70 kg	10 pieces

# Medical surface disinfection



## Beta Guard rfu

For the disinfection of Gel-a-Med® products



Beta Guard rfu organic wipes



Beta Guard rfu solution

Article-No.	Name	Description
M7-80502000DW	Beta Guard rfu organic wipes	Ready to use disinfectant wipes. Package contents: 80 pcs
M7-1000ML-OD	Beta Guard rfu solution	Ready-to-use spray/wipe disinfectant solution. Bottle contents: 1 liter



Beta Guard rfu products are CE marked and free of alcohol and aldehyde. They are used for the disinfection of surfaces of non-invasive medical products and surfaces in the medical environment. Dried surfaces do not stick after disinfection. Beta Guard products also show high material compatibility and are ideally suited for the disinfection of our Gel-a-Med® medical products.



## Scar plaster



### Gel-a-Med® Scar plaster

#### Silicone gel pads

Transparent and reusable. Our 100% medical grade silicone gel-based Gel-a-Med® scar patches are used to heal scar tissue and for the treatment and prevention of hypertrophic scars and keloids. The patch minimizes skin movement and helps to maintain the skin's natural moisture, thereby moisturizing the scar and improving regeneration.



Scar plaster

Article-No.	Description	PU
M1-03560003N	Scar plaster silicone gel	20 pieces
M1-0510003N	Scar plaster silicone gel	10 pieces
M1-0520003N	Scar plaster silicone gel	10 pieces
M1-1030003N	Scar plaster silicone gel	2 pieces (pair)
M1-12150003N	Scar plaster silicone gel	1 pieces
M1-2030003N	Scar plaster silicone gel	2 pieces (pair)


## Scar plaster application




1  
Measure



2  
Cutting



3  
Remove transparent  
film



4  
Apply



Make sure that a distance of at least 1,0 cm is maintained between the scar and the edge of the plaster on each side to ensure adequate coverage.

### Wearing time

- ✓ The optimum wearing time is 24 hours a day
- ✓ The minimum wearing time is 8-12 hours a day
- ✓ The recommended duration of treatment is at least one month

### Care

- ✓ Wash the plaster after each use with warm water and mild soap and allow it to air dry
- ✓ With the right care, you can extend the life of the plaster by 2 to 3 weeks

### Ingredients

- ✓ 100 % medical grade silicone
- ✓ Latex-free (latex-free plasters can therefore minimize many allergic reactions and have a positive effect on wound healing)

# Our service

## Free samples

Order your free Gel-a-Med® product sample now and experience the uniquely efficient pressure relief through silicone gel cushions. The sample is 150 x 100 mm in size and can be used, for example, as an ergonomic wrist cooling or heating pad.

Order free Gel-a-Med sample – here's how:

- ✓ Online at [www.gel-a-med.com](http://www.gel-a-med.com)
- ✓ By e-mail to [info@gel-a-med.com](mailto:info@gel-a-med.com)

### Your free sample

Made of silicone gel  
150 x 100 mm  
usable as an ergonomic wrist pad



## Special sizes · Special designs

Special sizes can be ordered according to individual requirements. With many years of experience in medical device manufacturing, we develop the perfect product solution together with our customers. We are happy to manufacture gel components according to your specifications for further processing or for sale under your own brand name.

The possibilities are virtually unlimited:

- ✓ Freely selectable dimensions
- ✓ Freely selectable height/thickness
- ✓ Individual shape
- ✓ Individual coloring or clear/transparent
- ✓ On request with cover linen (fabric / color freely selectable)



We will be happy to advise you and look forward to hearing from you!

**+49 211 556333 · [info@gel-a-med.com](mailto:info@gel-a-med.com)**

Also ask for our products in the fields of  
Living & Care · Automotive · Motorcycle · Aviation · Railway





# Gel-a-Med<sup>®</sup>

Medical sitting and lying

